

SMARCA2 / BRM Rabbit pAb

ABZ1600

Product Overview

Name	SMARCA2 / BRM Rabbit pAb
Accession(Primary)	P51531
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 181 kDa; Observed MW: 181 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P51531)

Gene Name SMARCA2

About Gene Symbol
(kD)



