

SOX17 Rabbit pAb

ABZ1601


Product Overview

Name	SOX17 Rabbit pAb
Accession(Primary)	Q9H6I2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 44 kDa; Observed MW: 44 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q9H6I2)

Gene Name SOX17

About Gene Symbol



Western blot analysis of SOX17 in Jurkat lysates using Sox17 antibody.

