

APOL6 Rabbit pAb

ABZ1602


Product Overview

Name	APOL6 Rabbit pAb
Accession(Primary)	Q9BWW8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q9BWW8)

Gene Name APOL6

About Gene Symbol



Western blot analysis of ApoL6 in K562 lysates using ApoL6 antibody.

