

NUP88 Rabbit pAb

ABZ1603


Product Overview

Name	NUP88 Rabbit pAb
Accession(Primary)	Q99567
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 84 kDa; Observed MW: 84 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q99567)

Gene Name NUP88

About Gene Symbol



Western blot analysis of NUP88 in K562 lysates using Nup88 antibody.

