

# Cytochrome P450 2A6/7/13 Rabbit pAb

ABZ1604


## Product Overview

|                             |                                            |
|-----------------------------|--------------------------------------------|
| <b>Name</b>                 | Cytochrome P450 2A6/7/13 Rabbit pAb        |
| <b>Accession(Primary)</b>   | P11509                                     |
| <b>Host Species</b>         | Rabbit                                     |
| <b>Clonality</b>            | Polyclonal                                 |
| <b>Isotype</b>              | IgG                                        |
| <b>Species Reactivity</b>   | Human                                      |
| <b>Calculated MW</b>        | Calculated MW: 56 kDa; Observed MW: 56 kDa |
| <b>Recommended Dilution</b> | WB-1:500-1:2000; ELISA-1:40000             |

## Target information (P11509)

**Gene Name** CYP2A6/CYP2A7/CYP2A13

**About Gene Symbol**



Western blot analysis of Cytochrome P450 2A6/7/13 in Jurkat lysates using CYP2A6/7/13 antibody.

