

NTT4 Rabbit pAb

ABZ1605

Product Overview

Name	NTT4 Rabbit pAb
Accession(Primary)	Q9H1V8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 81 kDa; Observed MW: 81 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000

Target information (Q9H1V8)

Gene Name SLC6A17**About Gene Symbol**

