

SCML2 Rabbit pAb

ABZ1607


Product Overview

Name	SCML2 Rabbit pAb
Accession(Primary)	Q9UQR0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 77 kDa; Observed MW: 77 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q9UQR0)

Gene Name SCML2

About Gene Symbol



Western blot analysis of SCML2 in 293 lysates using SCML2 antibody.

