

Wnt16 Rabbit pAb

ABZ1608


Product Overview

Name	Wnt16 Rabbit pAb
Accession(Primary)	Q9UBV4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:40000

Target information (Q9UBV4)

Gene Name WNT16

About Gene Symbol



Western blot analysis of Wnt16 in rat kidney lysates using Wnt16 antibody.

