

HspA6/A7 Rabbit pAb

ABZ1609

Product Overview

Name	HspA6/A7 Rabbit pAb
Accession(Primary)	P48741
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 40 kDa; Observed MW: 40 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P48741)

Gene Name HSPA7**About Gene Symbol**

Western blot analysis of HspA6/A7 in K562 lysates using HspA6/A7 antibody.

