

Synaptonemal Complex Protein 2 Rabbit pAb


ABZ1610

Product Overview

Name	Synaptonemal Complex Protein 2 Rabbit pAb
Accession(Primary)	Q9BX26
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 176 kDa; Observed MW: 176 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300; ELISA 1:2000-20000

Target information (Q9BX26)

Gene Name	SYCP2
About Gene Symbol	



Western blot analysis of Synaptonemal Complex Protein 2 in SKOV-3, BT474, Jurkat, PC-3, U937, MCF-7 lysates using Synaptonemal Complex Protein 2 antibody.