

TXNDC17 Rabbit pAb

ABZ1611


Product Overview

Name	TXNDC17 Rabbit pAb
Accession(Primary)	Q9BRA2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q9BRA2)

Gene Name TXNDC17

About Gene Symbol



Western blot analysis of Thioredoxin Domain Containing Protein 17 in HL-60, Raji, K562, SW480, A549, 22RV1 lysates using Thioredoxin Domain Containing Protein 17 antibody.