

Endothelin 1 Rabbit pAb

ABZ1612

Product Overview

Name	Endothelin 1 Rabbit pAb
Accession(Primary)	P05305
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P05305)

Gene Name	EDN1
About Gene Symbol	



Western blot analysis of Endothelin 1 in mouse lung lysates using Endothelin 1 antibody.