

HGS Rabbit pAb


ABZ1613

Product Overview

Name	HGS Rabbit pAb
Accession(Primary)	O14964
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 86 kDa; Observed MW: 86 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:20000

Target information (O14964)

Gene Name	HGS
About Gene Symbol	



Western blot analysis of HGS in mouse eye, mouse kidney, mouse testis lysates using HGS antibody.