


CD1c Rabbit pAb

ABZ1614

Product Overview

Name	CD1c Rabbit pAb
Accession(Primary)	P29017
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P29017)

Gene Name CD1C**About Gene Symbol**

Western blot analysis of CD1C in K562 lysates using CD1C antibody.

