

PSPN Rabbit pAb

ABZ1615

Product Overview

Name	PSPN Rabbit pAb
Accession(Primary)	O60542
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (O60542)

Gene Name PSPN**About Gene Symbol**

Western blot analysis of PSPN in 293 lysates using PSPN antibody.

