

Interferon epsilon Rabbit pAb

ABZ1616

Product Overview

Name	Interferon epsilon Rabbit pAb
Accession(Primary)	Q86WN2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 24 kDa; Observed MW: 24 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q86WN2)

Gene Name IFNE

About Gene Symbol



