

KIR2DL2/3 Rabbit pAb

ABZ1618

Product Overview

Name	KIR2DL2/3 Rabbit pAb
Accession(Primary)	P43628
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 38 kDa; Observed MW: 38 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P43628)

Gene Name KIR2DL3/KIR2DS2**About Gene Symbol**

Western blot analysis of KIR2DL2/3 in HeLa lysates using KIR2DL2/3 antibody.

