

Lymphotactin Rabbit pAb


ABZ1619

Product Overview

Name	Lymphotactin Rabbit pAb
Accession(Primary)	P47992
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 13 kDa; Observed MW: 13 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:20000

Target information (P47992)

Gene Name	XCL1/XCL2
About Gene Symbol	



Western blot analysis of Lymphotactin in K562
lysates using Lymphotactin antibody.