

CD85f Rabbit pAb

ABZ1620

Product Overview

Name	CD85f Rabbit pAb
Accession(Primary)	A6NI73
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (A6NI73)

Gene Name LILRA5**About Gene Symbol**

Western blot analysis of CD85f in HeLa lysates using CD85f antibody.

