


DCL-1 Rabbit pAb

ABZ1621

Product Overview

Name	DCL-1 Rabbit pAb
Accession(Primary)	Q8IX05
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q8IX05)

Gene Name CD302**About Gene Symbol**

Western blot analysis of DCL-1 in HeLa lysates using DCL-1 antibody.

