

Elongation Factor 1A1/2/P5 Rabbit pAb

ABZ1622


Product Overview

Name	Elongation Factor 1A1/2/P5 Rabbit pAb
Accession(Primary)	P68104
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 50 kDa; Observed MW: 50 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P68104)

Gene Name EEF1A1/EEF1A2/EEF1A1P5

About Gene Symbol



Western blot analysis of Elongation Factor 1A1/2/P5 in HepG2 lysates using Elongation Factor 1A1/2/P5 antibody.

