

IDO-2 Rabbit pAb

ABZ1623

Product Overview

Name	IDO-2 Rabbit pAb
Accession(Primary)	Q6ZQW0
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 47 kDa; Observed MW: 47 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (Q6ZQW0)

Gene Name IDO2

About Gene Symbol



