

Flt3 Ligand Rabbit pAb

ABZ1624

Product Overview

Name	Flt3 Ligand Rabbit pAb
Accession(Primary)	P49771
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 26 kDa; Observed MW: 26 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:20000

Target information (P49771)

Gene Name FLT3LG**About Gene Symbol**

