

CHRNA1 alpha1 Rabbit pAb

ABZ1625


Product Overview

Name	CHRNA1 alpha1 Rabbit pAb
Accession(Primary)	P02708
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 52 kDa; Observed MW: 55 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:5000; IF-1:50-1:200

Target information (P02708)

Gene Name CHRNA1

About Gene Symbol



Western blot analysis of AChR alpha1 in NIH3T3 lysates using AChR alpha1 antibody.

