

STRA8 Rabbit pAb


ABZ1626

Product Overview

Name	STRA8 Rabbit pAb
Accession(Primary)	Q7Z7C7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 37 kDa; Observed MW: 66 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q7Z7C7)

Gene Name	STRA8
About Gene Symbol	



Western blot analysis of STRA8 in NIH3T3, 293 lysates using STRA8 antibody.

