

NOXIN Rabbit pAb

ABZ1627


Product Overview

Name	NOXIN Rabbit pAb
Accession(Primary)	Q8IXT1
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 112 kDa; Observed MW: 85 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q8IXT1)

Gene Name DDIAS

About Gene Symbol



Western blot analysis of NOXIN in 823, 293T, AD293, Hela, HepG2 lysates using NOXIN antibody.

