

Tensin 2 Rabbit pAb

ABZ1628

Product Overview

Name	Tensin 2 Rabbit pAb
Accession(Primary)	Q63HR2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 153 kDa; Observed MW: 125 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q63HR2)

Gene Name TNS2**About Gene Symbol**

