


ASAH3 Rabbit pAb

ABZ1629

Product Overview

Name	ASAH3 Rabbit pAb
Accession(Primary)	Q8TDN7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 31 kDa; Observed MW: 34 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q8TDN7)

Gene Name ACER1**About Gene Symbol**

Western blot analysis of ASAH3 in HepG2 lysates using ASAH3 antibody.

