


TCEAL6 Rabbit pAb

ABZ1630

Product Overview

Name	TCEAL6 Rabbit pAb
Accession(Primary)	Q6IPX3
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 22 kDa; Observed MW: 32 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q6IPX3)

Gene Name TCEAL6**About Gene Symbol**

Western blot analysis of TCEAL6 in CoLo lysates using TCEAL6 antibody.

