


ANKRD20A1 Rabbit pAb

ABZ1631

Product Overview

Name	ANKRD20A1 Rabbit pAb
Accession(Primary)	Q5TYW2
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 94 kDa; Observed MW: 94 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q5TYW2)

Gene Name ANKRD20A1**About Gene Symbol**

Western blot analysis of ANKRD20A1 in 3T3
lysates using ANKRD20A1 antibody.

