

TRMT11 Rabbit pAb

ABZ1632


Product Overview

Name	TRMT11 Rabbit pAb
Accession(Primary)	Q7Z4G4
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 53 kDa; Observed MW: 53 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q7Z4G4)

Gene Name TRMT11

About Gene Symbol



Western blot analysis of TRMT11 in HepG2 , JK, Colo, Hela lysates using TRMT11 antibody.

