

PIK3C2G Rabbit pAb

ABZ1633


Product Overview

Name	PIK3C2G Rabbit pAb
Accession(Primary)	O75747
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 166 kDa; Observed MW: 160 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (O75747)

Gene Name PIK3C2G

About Gene Symbol



Western blot analysis of PI3 Kinase C2 gamma in A549 lysates using PI3 Kinase C2 gamma antibody.

