


PODXL2 Rabbit pAb

ABZ1634

Product Overview

Name	PODXL2 Rabbit pAb
Accession(Primary)	Q9NZ53
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 65 kDa; Observed MW: 65 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (Q9NZ53)

Gene Name PODXL2**About Gene Symbol**

Western blot analysis of Endoglycan in Hela lysates using Endoglycan antibody.

