

Cripto 3 Rabbit pAb

ABZ1635


Product Overview

Name	Cripto 3 Rabbit pAb
Accession(Primary)	P51864
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 21 kDa; Observed MW: 21 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000

Target information (P51864)

Gene Name TDGF1P3

About Gene Symbol



Western blot analysis of Cripto 3 in HepG2, 293T ,
AD293 , HeLa lysates using Cripto 3 antibody.

