

Collagen XXVIII alpha 1 Rabbit pAb

ABZ1636

Product Overview

Name	Collagen XXVIII alpha 1 Rabbit pAb
Accession(Primary)	Q2UY09
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 117 kDa; Observed MW: 117 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q2UY09)

Gene Name COL28A1**About Gene Symbol**

