

Plasminogen Receptor Rabbit pAb

ABZ1637

Product Overview

Name	Plasminogen Receptor Rabbit pAb
Accession(Primary)	Q9HBL7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Calculated MW	Calculated MW: 17 kDa; Observed MW: 17 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q9HBL7)

Gene Name PLGRKT**About Gene Symbol**

