

Placental Lactogen Rabbit pAb


ABZ1638

Product Overview

Name	Placental Lactogen Rabbit pAb
Accession(Primary)	P01243
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 25 kDa; Observed MW: 25 kDa
Recommended Dilution	WB-1:500-1:2000; IHC:-1:100-1:300; IF-1:50-1:200; ELISA-1:10000-1:20000

Target information (P01243)

Gene Name	CSH1
About Gene Symbol	



Western blot analysis of Placental Lactogen in mouse brain lysates using Placental Lactogen antibody.