


Cactin Rabbit pAb

ABZ1639

Product Overview

Name	Cactin Rabbit pAb
Accession(Primary)	Q8WUQ7
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 89 kDa; Observed MW: 95 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q8WUQ7)

Gene Name CACTIN**About Gene Symbol**

Western blot analysis of Cactin in hepg2 lysates using Cactin antibody.

