

HLA A Rabbit pAb


ABZ1640

Product Overview

Name	HLA A Rabbit pAb
Accession(Primary)	P04439
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Calculated MW	Calculated MW: 41 kDa; Observed MW: 41 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (P04439)

Gene Name	HLA-A
About Gene Symbol	HLA-A (Major Histocompatibility Complex, Class I, A) is a membrane receptor involved in adaptive immunity. It is an important research target in cell signaling, drug discovery, and cancer research and is widely studied to elucidate cellular mechanisms and



Western blot analysis of HLA A in HELA lysates
using HLA A antibody.