

Cytokeratin Rabbit pAb

ABZ1641


Product Overview

Name	Cytokeratin Rabbit pAb
Accession(Primary)	P35908
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 65 kDa; Observed MW: 70 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P35908)

Gene Name

Keratin

About Gene Symbol


Western blot analysis of Cytokeratin in mouse lung, mouse brain, mouse kidney, mouse heart lysates using Cytokeratin antibody.

