


RPA14/RPA3 Rabbit pAb

ABZ1642

Product Overview

Name	RPA14/RPA3 Rabbit pAb
Accession(Primary)	P35244
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 14 kDa; Observed MW: 14 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P35244)

Gene Name RPA3**About Gene Symbol**

Western blot analysis of RPA3 in SW480, MCF-7 lysates using RPA3 antibody.

