


Macrophage Inflammatory Protein 5 Rabbit pAb

ABZ1644

Product Overview

Name	Macrophage Inflammatory Protein 5 Rabbit pAb
Accession(Primary)	Q16663
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (Q16663)

Gene Name CCL15**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human liver using Macrophage Inflammatory Protein 5 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.