

GRO gamma Rabbit pAb

ABZ1645


Product Overview

Name	GRO gamma Rabbit pAb
Accession(Primary)	P19876
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:100-1:300; ELISA-1:10000; IF-1:50-1:200

Target information (P19876)

Gene Name CXCL3

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human heart using GRO gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.