


IL-17B Rabbit pAb

ABZ1646

Product Overview

Name	IL-17B Rabbit pAb
Accession(Primary)	Q9UHF5
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (Q9UHF5)

Gene Name IL17B**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human small intestine using IL-17B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.