


AVPR1B Rabbit pAb

ABZ1647

Product Overview

Name	AVPR1B Rabbit pAb
Accession(Primary)	P47901
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (P47901)

Gene Name AVPR1B**About Gene Symbol**

Immunohistochemistry analysis of paraffin-embedded Human brain using AVPR1B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.