

Insulin Rabbit pAb

ABZ1648


Product Overview

Name	Insulin Rabbit pAb
Accession(Primary)	P01308
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Recommended Dilution	IHC-1:50-1:200; ELISA-1:10000-1:20000; IF-1:50-1:200

Target information (P01308)

Gene Name INS

About Gene Symbol



Immunohistochemistry analysis of paraffin-embedded Human pancreas using Insulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.