


NUT Rabbit pAb

ABZ1649

Product Overview

Name	NUT Rabbit pAb
Accession(Primary)	Q86Y26
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 120 kDa; Observed MW: 63 kDa
Recommended Dilution	WB-1:500-1:2000; ELISA-1:10000-1:20000

Target information (Q86Y26)

Gene Name NUTM1**About Gene Symbol**

Western blot analysis of NUT in 293T lysates
using NUT antibody.

