


KCTD7 Rabbit pAb

ABZ1651

Product Overview

Name	KCTD7 Rabbit pAb
Accession(Primary)	Q96MP8
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse, Rat
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (Q96MP8)

Gene Name KCTD7**About Gene Symbol**

Western blot analysis of KCTD7 in mouse brain ,
mouse spleen lysates using KCTD7 antibody.

