


Integrin beta 7 Rabbit pAb

ABZ1652

Product Overview

Name	Integrin beta 7 Rabbit pAb
Accession(Primary)	P26010
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 87 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:50-1:300

Target information (P26010)

Gene Name ITGB7**About Gene Symbol**

Western blot analysis of Integrin beta 7 in mouse heart lysates using Integrin beta 7 antibody.

