

Collagen XIX alpha 1 Rabbit pAb


ABZ1653

Product Overview

Name	Collagen XIX alpha 1 Rabbit pAb
Accession(Primary)	Q14993
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Species Reactivity	Human, Mouse
Calculated MW	Calculated MW: 115 kDa; Observed MW: 130 kDa
Recommended Dilution	WB-1:500-1:2000; IHC-1:100-1:300; IF-1:200-1:1000; ELISA-1:10000

Target information (Q14993)

Gene Name	COL19A1
About Gene Symbol	



Western blot analysis of Collagen XIX alpha 1 in
HELA lysates using Collagen XIX alpha 1 antibody.